

AP4 Retrieve Release Notes

AP4 Retrieve 1.0

19 October 2006

Contents

ABOUT AP4 RETRIEVE	2
INSTALLATION	2
<i>System Requirements</i>	<i>2</i>
<i>To Install AP4 Retrieve</i>	<i>2</i>
<i>To Start AP4 Retrieve</i>	<i>2</i>
REVISION HISTORY	2
<i>AP4 Retrieve 1.0</i>	<i>2</i>

About AP4 Retrieve

AP4 Retrieve is the software for downloading data from the AP4 Porometer.

Installation

System Requirements

To install and use AP4Retrieve you will need:

- A PC running Windows 98, 2000, XP, or a later version.
- One free RS232 serial port, or USB-RS232 adapter.
- 10MB of available hard disk space.
- AP4 RS232 cable (supplied with the AP4).

To Install AP4 Retrieve

You need local administrator rights to install AP4 Retrieve 1.0 on computers with Windows 2000, XP and Vista. Please contact your system administrator for further information.

Select the **Install AP4 Retrieve 1.0** link in the Contents page of the Delta-T Software and Manuals CD and follow the on-screen instructions.

Alternatively, download AP4 Retrieve 1.0 from www.delta-t.co.uk and click **Run** when the download completes.

The installation program creates an **AP4 Retrieve** shortcut in the Windows Programs menu.

To Start AP4 Retrieve

Select the **AP4 Retrieve** shortcut in the Windows Programs menu.

Revision History

AP4 Retrieve 1.0

AP4Retrieve 1.0 is the first release of AP4Retrieve.